



WHY TOP PROFESSIONALS CHOOSE CTi KNEE BRACES FOR MX AND FMX

CTi braces are medical grade products, covered by most insurance plans with a doctor's prescription. Whether custom or off-the-shelf, every CTi ligament brace benefits from the input of medical professionals and feedback from professional riders. The resulting brace is the perfect blend of support and function, designed to help you rehab and ride safely after a knee injury.

MX and FMX riders have unique demands and anatomy. CTi features specific design advantages and/or adjustments that can be made to optimize performance for riding dirt bikes, including:

- A basic frame design that functions well for riding. The CTi is the only ligament knee brace that captures the tibia, locking down securely on the tibial crest. The rigid upper and lower arms combine with the hinges, condyle pads and straps to properly align the knee joint and provide an accurately-tracking exoskeletal support.
- There are no rigid components in the upper medial aspect. This keeps the rider streamlined to the bike and enables them to use their thighs to better grip the tank without pinching or digging in.
- CTi can be ordered with a patella cup. This protective cup keeps the patella safe from flying rocks and other debris common to riding environments.
- CTi has no rigid components across the posterior (back) side of the brace, so you won't experience pinching or cramping when in deep flexion (squatting) through whoops or tight corners .
- CTi's low-profile design allows it to tuck into the top of MX boots. It actually works in conjunction with the boot to help provide better overall protection.
- CTi can be manufactured to be shorter, longer, thicker or thinner based upon your height and weight.

FOR MORE INFORMATION AND HOW TO ORDER, CALL 800.233.6263 OR VISIT OSSUR.COM/CTi