

ÖSSUR®

FORM FIT®

THE CHOICE OF EXPERTS



Form Fit Ankle and Form Fit Ankle with Figure 8

Sometimes, our bodies need assistance. For decades, specialists have trusted Form Fit supports, a complete range of quality products by Össur. A global leader in innovative orthopaedic products, Össur's biomechanical expertise has helped millions live Life Without Limitations.

That's why Form Fit is the Choice of Experts.



FORMFIT®

THE CHOICE OF EXPERTS



FORM FIT ANKLE AND FORM FIT ANKLE WITH FIGURE 8

The Form Fit Ankle brace provides stability and protection for the ankle – whether used prophylactically, for chronic instability or following an injury. A hybrid between a rigid ankle stirrup and a soft ankle support, the Form Fit Ankle ensures a close fit which limits inversion, eversion, flexion and extension movements. The low profile design fits easily into shoes and the quick lace up feature makes it easy to don and doff. Breathable material ensures it can be comfortably worn all day.

The Form Fit Ankle with Figure 8 strap version, offers an even greater level of support, by further limiting the range of motion of the ankle.

Features and Benefits

- Breathable material
- Quick lace Up
- Inversion/Eversion Control
- Figure 8 strap version

Indications

- Grade I-II sprains and strains
- Chronic instability
- Prophylactic use

ORDERING INFORMATION

FORM FITANKLE

Standard Part#	Figure 8 Part#	Size	Ankle Circumference
B-212010001	B-212000001	XSmall	11-12" (28-30.5cm)
B-212010002	B-212000002	Small	12-13" (30.5-33cm)
B-212010003	B-212000003	Medium	13-14" (33-35.5cm)
B-212010004	B-212000004	Large	14-15" (35.5-38cm)
B-212010005	B-212000005	XLarge	15-16" (38-40.5cm)

* Measure ankle circumference diagonally over the base of the heel and across the ankle bone.

FOLLOW ÖSSUR ON



TEL (800) 233-6263
CANADA (800) 663-5982
WWW.OSSUR.COM

