



FORM FIT® KNEE MCL

The Form Fit Knee MCL provides external stability and protection for low to medium contact/impact activities and while recovering from injury. A unique, single hinged, hybrid wrap/sleeve design sets this brace apart from the competition. The lateral, range of motion hinge keeps the medial side free of obstruction. The Dynamic Force Strap aids in the resistance of valgus forces during knee extension, while the adjustable tibial strap improves suspension by minimizing slippage.

- Single, steel ROM lateral hinge (0° or -20° extension)
- Dynamic Force Strap resist valgus forces
- Hybrid wrap/sleeve for easy application
- Adjustable tibial strap for improved suspension

INDICATIONS

- Mild to Moderate Medial Collateral Ligament (MCL) sprains (Grade I and Grade II)

PART#	SIZE	THIGH CIRCUMFERENCE
B-237X00001	XSmall	13-15.5" (33-39cm)
B-237X00002	Small	15.5-17.5" (39-44cm)
B-237X00003	Medium	17.5-20.5" (44-52cm)
B-237X00004	Large	20.5-23.5" (52-60cm)
B-237X00005	XLarge	23.5-26.5" (60-67cm)

X = 5 (Right) or 6 (Left)

* Measure thigh circumference 6" above mid-patella