

ÖSSUR®

FORMFIT®

THE CHOICE OF EXPERTS



Form Fit Shoulder Brace with Abduction

Sometimes, our bodies need assistance. For decades, specialists have trusted Form Fit supports, a complete range of quality products by Össur. A global leader in innovative orthopaedic products, Össur's biomechanical expertise has helped millions live Life Without Limitations.

That's why Form Fit is the Choice of Experts.



FORM FIT®

THE CHOICE OF EXPERTS



FORM FIT SHOULDER BRACE WITH ABDUCTION

The Form Fit® Shoulder Brace is designed to provide immobilization following shoulder injury/surgery. With its streamlined fit, sporty look and breathability it offers excellent comfort and ease of use.

Same as its predecessor, the SmartSling®, the Form Fit Shoulder Brace also provides abduction of the shoulder (when the abduction pillow is attached), or external rotation of the shoulder (when the external rotation pillow is attached).

Indications

- Rotator cuff repairs
- Bankart repairs
- Capsular shifts
- Glenohumeral dislocations/subluxations
- Soft tissue repairs/strains
- Shoulder instabilities/injuries
- Post-surgical rehabilitation

Product Highlights

- Breathable material
- Quick-release buckles
- Offloading strap
- Extra padding
- Universal thumb loop
- 30° external rotation (when external rotation pillow is attached)

Optional axilla pillow helps distribute weight over a larger surface area and decreases pressure on the ulnar nerve.

ORDERING INFORMATION

FORM FIT SHOULDER BRACE WITH ABDUCTION

Part#	Size	Length
SSS	Small	10.5-13" (up to 29cm)
SSR	Medium	13.25-15.75" (29-40cm)
SSL	Large	16-17.5" (40-44.5cm)

Olecranon to end of closed fist

FOLLOW ÖSSUR ON



USA
(800)233-6263
WWW.OSSUR.COM

CANADA
(800)663-5982
WWW.OSSUR.CA

